



Minisart® RC Syringe Filters for HPLC Sample Preparation


RC versus PVDF used with Methanol and Acetonitrile

Comparison with PVDF


PVDF is a membrane material often used for HPLC sample preparation. The fluoro-chemical process needed to produce this raw material can be harmful to the environment. Moreover, PVDF requires an additional coating step with methacrylates to make its membrane surface hydrophilic. In addition, unlike RC syringe filters, those with PVDF membranes are not suitable for DMSO and other amides, or for ketones, esters, and ethers.

Conclusion

Minisart® RC is just as suitable, and often even better, for HPLC sample preparation than PVDF. The broader chemical compatibility of Minisart® RC enables usage with DMSO and for an even wider range of applications.



- Methanol filtered using PVDF, Competitor B
- Methanol filtered using PVDF, Competitor A
- Methanol filtered using Minisart® RC
- Unfiltered Methanol



- Acetonitrile filtered using PVDF, Competitor B
- Acetonitrile filtered using PVDF, Competitor A
- Acetonitrile filtered using Minisart® RC
- Unfiltered Acetonitrile

HPLC Procedure

Column C18: 250 × 4.6 mm, Flow Rate: 1 ml/min, Wavelength: 210 nm

HPLC Injection Volume: 20 µl, Analysis Time: 65 min, Temperature: 40°C, Mobile Phases: A) Acetonitrile | B) Water, Gradient: Hold 60% A for 10 min | 60% to 95% A in 20 min | 95% to 100% A in 35 min

Reliable and Ultrapure

Clean & Green

Are you using high-purity filter paper in your lab? Our regenerated cellulose (RC) membrane incorporated in Minisart® RC syringe filters is just as pure, but all-purpose and more convenient.

RC is produced using renewable raw materials in an eco-friendly process, which uses biodegradable components and recycles all solvents employed during manufacture.

Our RC polymer is naturally hydrophilic so no additional coating step is necessary. Try our Clean and Green solution to prepare your samples.



Ordering Information

Ø mm	Pore Size	Sterile*	Qty/Pk	Order No.
Minisart® RC (Regenerated Cellulose + PP)				
25 mm	0.2 µm	Yes	50	17764-----ACK
25 mm	0.2 µm	No	50	17764-----K
25 mm	0.2 µm	No	200	17764-----S
25 mm	0.2 µm	No	500	17764-----Q
25 mm	0.45 µm	No	50	17765-----K
25 mm	0.45 µm	No	200	17765-----S
25 mm	0.45 µm	No	500	17765-----Q
15 mm	0.2 µm	Yes	50	17761-----ACK
15 mm	0.2 µm	No	50	17761-----K
15 mm	0.2 µm	No	500	17761-----Q
15 mm	0.45 µm	No	50	17762-----K
15 mm	0.45 µm	No	500	17762-----Q
4 mm	0.2 µm	No	50	17821-----K
4 mm	0.2 µm	No	500	17821-----Q
4 mm	0.45 µm	No	50	17822-----K
4 mm	0.45 µm	No	500	17822-----Q

* Sterile Minisart® are individually packaged and have been sterilized by ethylene oxide (EO).
Not presterilized Minisarts can be sterilized by autoclaving at 121°C for 30 min | or by using EO.

Would You Like to Use Other Membranes?

Please refer to Minisart® NY, or Minisart® SRP for highest chemical compatibility.

Do You Need Minisart® with Pre-filters for Filtration of Highly Particle Laden Samples?

Please refer to Minisart® NY Plus which incorporate a glass fiber pre-filter.

Are You Looking for CE-certified Minisart®?

Please request CE-certified article numbers for Minisart® NML, Minisart® HY and Minisart® SRP.

Sartorius offers a full range of syringe filters dedicated for various filtration applications.

Sartorius Stedim Biotech GmbH
August-Spindler-Strasse 11
37079 Goettingen, Germany

Phone +49.551.308.0
Fax +49.551.308.3289



▶ www.sartorius.com

Specifications subject to change without notice.
Printed in Germany on paper that has been bleached without any use of chlorine. | W
Publication No.: SL-0016-e141001
Order No.: 85037-547-90 · Ver. 10 | 2014