Operating Instructions | Betriebsanleitung | Mode d'emploi | Instrucciones de manejo |  
Manuale d'uso | Руководство по эксплуатации | 操作说明 | 取扱説明書

Original Operating Instructions | Original-Betriebsanleitung | Mode d'emploi original |  
Instrucciones de manejo originales | Manuale d'uso originale | Оригинальное руководство  
по эксплуатации | 原始操作说明 | 取扱説明書 (オリジナル版)

# STAT-PEN

YSTP01

Device for neutralizing electrostatic charge | Gerät zur Neutralisierung elektrostatischer Aufladung |  
Appareil pour la neutralisation des charges électrostatiques | Aparato para neutralizar la carga electrostática |  
Apparecchio per neutralizzare le cariche elettrostatiche | Устройство для нейтрализации  
электростатического заряда | 静电荷中和设备 | 静电荷中和用機器

1000101778

SARTORIUS

English	page	3
Deutsch	Seite	16
Français	page	29
Español	página	42
Italiano	pagina	55
Русский	страница	68
中文	页	81
日本語	ページ	94

# Contents

<b>1</b>	<b>About these Instructions</b>	<b>4</b>	<b>7</b>	<b>Cleaning and Maintenance</b>	<b>10</b>
1.1	Validity	4	7.1	Cleaning the Device	10
1.2	Target Groups	4	7.2	Maintenance Schedule	10
1.3	Symbols Used	4	7.3	Restarting	10
1.3.1	Warnings in Operation Descriptions	4	<b>8</b>	<b>Malfunctions</b>	<b>11</b>
1.3.2	Other Symbols	4	8.1	LED	11
<b>2</b>	<b>Safety Information</b>	<b>5</b>	8.2	Checking the Operation of the Device	11
2.1	Intended Use	5	<b>9</b>	<b>Decommissioning</b>	<b>12</b>
2.1.1	Modifications to the Device	5	9.1	Decommission the Device	12
2.1.2	Repairs to the Device	5	<b>10</b>	<b>Storage and Shipping</b>	<b>12</b>
2.2	Personnel Qualification	5	10.1	Decontaminating the Device	12
2.3	Significance of these Instructions	5	10.2	Storage	12
2.4	Device Functionality	5	10.3	Returning the Device and Parts	12
2.5	Electrical Equipment	6	<b>11</b>	<b>Disposal</b>	<b>12</b>
2.5.1	Power Supply Cable	6	11.1	Disposing of the Device and Parts	12
2.6	Conduct in an Emergency	6	<b>12</b>	<b>Technical Data</b>	<b>13</b>
2.7	Accessories, Consumables, and Spare Parts	6	12.1	Dimensions and Weights	13
2.8	Danger of Injury due to Defective Device Components	6	12.2	Ambient Conditions	13
2.8.1	Damage to the Housing	6	12.3	Electrical Specifications	13
2.8.2	Damage to the Power Supply Cable	6	12.3.1	Power Supply	13
2.8.3	Damage to the Housing Insulation	6	12.3.2	Safety and Electromagnetic Compatibility	14
2.9	Ingress of Fluids	6	12.3.3	Neutralizing	14
<b>3</b>	<b>Device Description</b>	<b>7</b>	12.4	Operating Area	14
3.1	Device Overview	7	12.5	Discharging Times	15
3.2	LEDs	7	12.6	Cleaning Agents and Cleaning Processes	15
3.3	Symbols on the Device	7	<b>13</b>	<b>Sartorius Service</b>	<b>15</b>
3.3.1	Warning Symbols	7	<b>14</b>	<b>Conformity Documents</b>	<b>15</b>
3.3.2	Warnings	7			
<b>4</b>	<b>Installation</b>	<b>8</b>			
4.1	Scope of Delivery	8			
4.2	Selecting an Installation Site	8			
4.3	Unpacking	8			
4.4	Acclimatizing	8			
<b>5</b>	<b>Getting Started</b>	<b>9</b>			
5.1	Connecting the Power Supply	9			
<b>6</b>	<b>Operation</b>	<b>9</b>			
6.1	Switching the Device On or Off	9			
6.2	Neutralizing Electrostatic Charge on Objects	9			

# 1 About these Instructions

## 1.1 Validity

These instructions are part of the device. These instructions apply to the following version of the device:

Device	Type
STAT-PEN	YSTP01

## 1.2 Target Groups

These instructions are designed for the following target groups. The target groups must possess the knowledge specified below.

Target group	Knowledge and Qualifications
User	The user is responsible for operating the device and for the associated work processes. The user understands the hazards which may arise when working with the device and knows how to prevent them. The user has been trained to operate the device.
Operator	The operator of the device is responsible for ensuring compliance with workplace health and safety regulations. The operator must ensure that all persons who work with the device have access to the relevant information and are trained in working with the device.

## 1.3 Symbols Used

### 1.3.1 Warnings in Operation Descriptions

#### CAUTION

Denotes a danger that can result in moderate or minor injury if it is **not** avoided.

### 1.3.2 Other Symbols

- ▶ Required action: Describes activities that must be carried out. The actions in the sequence must be carried out in succession.
- ▷ Result: Describes the result of the actions carried out.

## 2 Safety Information

### 2.1 Intended Use

The device is used to neutralize electrostatic charge on objects to which the ionization pin is pointed. The ionization pin generates positive and negative ions.

The device is intended exclusively for use in accordance with these instructions. Any other use is considered **improper**.

If the device and accessories are used **improperly**: The device's protective systems may be impaired. This can lead to unforeseeable personal injury.

#### Operating Conditions for the Device

Do **not** use the device in potentially explosive environments. Only use the device indoors.

The device may only be used with the equipment and under the operating conditions described in the Technical Data section of these instructions.

#### 2.1.1 Modifications to the Device

If the device is modified, for example by attaching extra components: The safety of the device may be impaired or the device compliance may lose its validity.

If you have any queries regarding modifications to the device, contact Sartorius.

#### 2.1.2 Repairs to the Device

Device repairs may only be carried out by persons with appropriate specialized knowledge of the device. If the device is **not** repaired by a specialist: The safety of the device may be impaired or the test marks may lose their validity.

We recommend that any repair work, even that not covered by the warranty, is carried out by Sartorius Service or after consulting with Sartorius Service.

### 2.2 Personnel Qualification

If individuals who do **not** have sufficient knowledge on the safe handling of the device carry out work on the device: Those individuals may injure themselves or other people nearby.

- ▶ Ensure that all persons working on the device possess the necessary knowledge and qualifications (for description, see chapter "1.2 Target Groups", page 4).
- ▶ If a particular qualification is required for the actions described: Have these activities carried out by the required target group.
- ▶ If **no** particular qualification is required for the actions described: Have these activities carried out by the "user" target group.

### 2.3 Significance of these Instructions

Failure to follow the instructions in this manual may have serious consequences, e.g. danger to individuals.

- ▶ Read the instructions carefully and completely. The required actions in the instructions build on each other.
- ▶ Ensure that the information contained in these instructions is available to all individuals working on the device.
- ▶ Retain the instructions.
- ▶ If these instructions are lost, request a replacement or download the latest version from the Sartorius website ([www.sartorius.com](http://www.sartorius.com)).

### 2.4 Device Functionality

A damaged device or worn-out parts can cause malfunctions or lead to hard-to-detect hazards.

- ▶ Only operate the device when it is safe and in proper working order.
- ▶ Have any malfunctions or damage repaired immediately by the Sartorius Service.

## 2.5 Electrical Equipment

### 2.5.1 Power Supply Cable

Serious injury can result, e.g. from electric shocks, if an unsuitable/inadequately dimensioned power supply cable is used.

- ▶ Only use the original power supply cable.
- ▶ If the power supply cable needs to be replaced: Contact Sartorius Service. Do **not** repair or modify the power supply cable.

## 2.6 Conduct in an Emergency

If an emergency occurs, e.g., due to malfunctions of the device or dangerous situations: People may be injured. The device must be immediately taken out of operation:

- ▶ Disconnect the device from the power supply.
- ▶ Secure the device to prevent it from restarting.

## 2.7 Accessories, Consumables, and Spare Parts

Unsuitable accessories, consumables, and spare parts can adversely affect functionality and safety and can have the following consequences:

- Risk of injury to persons
  - Damage to the device
  - Device malfunctions
  - Failure of the device
- ▶ Only use approved accessories, consumables, and spare parts supplied by Sartorius.
  - ▶ Only use accessories, consumables, and spare parts that are in technically perfect condition.

## 2.8 Danger of Injury due to Defective Device Components

### 2.8.1 Damage to the Housing

In case of contact with damaged housing during operation: Individuals may suffer an electric shock. This can lead to life-threatening injuries.

- ▶ If the housing is damaged: Disconnect the device from the power supply and dispose of it.

### 2.8.2 Damage to the Power Supply Cable

If a damaged power supply cable is used or in case of contact with a damaged power supply cable: Individuals may suffer an electric shock. This can lead to life-threatening injuries.

- ▶ If the power supply cable is damaged: Replace the power supply cable. For a replacement, contact Sartorius Service. Do **not** repair or modify the power supply cable.

### 2.8.3 Damage to the Housing Insulation

If the housing insulation is damaged due to high temperatures: Individuals may suffer an electric shock. This can lead to life-threatening injuries.

- ▶ Only operate the device at the permissible temperatures and altitudes (see chapter “12.2 Ambient Conditions”, page 13).

## 2.9 Ingress of Fluids

If fluids enter the device: Individuals may suffer an electric shock. This can lead to life-threatening injuries.

- ▶ If fluids have entered the device:
  - ▶ Disconnect the device from the power supply.
  - ▶ Secure the device to prevent it from restarting.

## 3 Device Description

### 3.1 Device Overview


Fig. 1: Front view of the device

Pos.	Name	Description
1	On   Off switch	<ul style="list-style-type: none"> <li>Position 1: Switch on the device.</li> <li>Position 0: Switch off the device.</li> </ul>
2	Ionization pin holder	
3	Power supply connection	
4	Ionization pin connection	The connection cable is securely connected.
5	Manufacturer's ID label	
6	Ionization microswitch	If held down: The ionization is activated.
7	Ionization pin	
8	Ion emitter	4 ion emitters generate positive and negative ions.

### 3.2 LEDs


Fig. 2: LEDs

Pos.	Name	Description
1	LED	No function.
2	Power supply LED	Lit red: The device is switched on.
3	Ionization LED	Lit red: The ionization pin is activated.

### 3.3 Symbols on the Device

#### 3.3.1 Warning Symbols

Symbol	Meaning
	Electrical voltage: Risk of injury due to live parts. Only electricians may have access to and work on these parts.
	Hazards on the device: Risk of injury due to device components
	Power supply: The device is connected to the power supply and switched on.

#### 3.3.2 Warnings

Warning	Meaning
<b>CAUTION</b> Risk of electric shock Do not open	<b>CAUTION</b> Risk of electrocution Do not open the device

## 4 Installation

### 4.1 Scope of Delivery

Item	Quantity
Device	1
Power supply cable, fixed to the device	1
Ionization pin	1
Operating instructions	1

### 4.2 Selecting an Installation Site

#### Procedure

- ▶ Make sure that the following conditions are met at the installation site:

Condition	Requirements
Ambient conditions	– Suitability tested (see chapter 12.2, page 13 for ambient conditions)
Setup surface	– Stable and level surface – Sufficiently dimensioned for the device and the operating area (for device space requirements see chapter “12.1 Dimensions and Weights”, page 13 and chapter “12.4 Operating Area”, page 14).
Access to operation-relevant parts	Convenient and safe

### 4.3 Unpacking

#### Procedure

- ▶ Unpack the device.
- ▶ **⚠ CAUTION** Risk of crushing if the device falls! Make sure that the device has a suitable setup surface (see chapter “4.2 Selecting an Installation Site”, page 8).
- ▶ Install the device at the intended installation site.

### 4.4 Acclimatizing

When a cold device is brought into a warm environment: the temperature difference can lead to condensation (moisture formation) from humidity in the device. Moisture in the device can lead to malfunctions.

#### Procedure

- ▶ Allow the device to acclimatize for approx. 2 hours at the installation site. The device must be disconnected from the power supply during that time.

## 5 Getting Started

### 5.1 Connecting the Power Supply

#### Procedure

- ▶ Check whether the voltage specifications on the manufacturer's ID label match those of the power supply at the installation site.
  - ▶ If the input voltage is too high: **Do not** connect the device to the power supply.
  - ▶ Contact Sartorius Service.
- ▶ Check whether the country-specific power plug matches the mains sockets at the installation site.
  - ▶ If necessary: Contact Sartorius Service.
- ▶ Set the On | Off switch to the 0 position.
- ▶ Connect the mains plug to the power supply at the installation site.

## 6 Operation

### 6.1 Switching the Device On or Off

#### Procedure

- ▶ To switch the device on: Set the On | Off switch to the 1 position.
  - ▷ The power supply LED lights up red.
- ▶ To switch the device off: Set the On | Off switch to the 0 position.

### 6.2 Neutralizing Electrostatic Charge on Objects

Positive and negative electrostatic charge can be neutralized.

#### Requirements

The device is switched on.

#### Procedure

- ▶ Remove the ionization pin from the ionization pin holder.
- ▶ Point the ionization pin at a suitable distance from the object to be neutralized (see chapter "12.4 Operating Area", page 14). The ion emitters must point to the object.
- ▶ Hold down the ionization microswitch.
  - ▷ The ionization pin is activated and both ionization LEDs light up.
  - ▷ The object's electrostatic charge is neutralized. The duration depends on the distance between the ionization pin and the object (see chapter "12.5 Discharging Times", page 15).
- ▶ To end the neutralization: Release the microswitch.
- ▶ Place the ionization pin in the ionization pin holder.

# 7 Cleaning and Maintenance

## 7.1 Cleaning the Device

### Procedure

- ▶ Switch off the device.
- ▶ Disconnect the device from the power supply at the installation site.
- ▶ Only use suitable cleaning agents and cleaning procedures and observe the product information for the cleaning agent used (see chapter "12.6 Cleaning Agents and Cleaning Processes", page 15).
- ▶ Protect the device against moisture and liquids.
- ▶ Wipe down the housing with a dry or slightly moistened fabric dust cloth.

## 7.2 Maintenance Schedule

Interval	Component	Activity	Chapter, page	Target group
Annually	Device	Have an electrician carry out a general inspection (in accordance with local government regulations).		Operator

## 7.3 Restarting

### Procedure

- ▶ Restart the device (see chapter "5 Getting Started", page 9).

## 8 Malfunctions

### 8.1 LED

Malfunction	Cause	Correction	Chapter, page
The power supply LED is <b>not</b> lighting up.	The device is switched off.	Switch on the device.	8.2, 11
	The device is defective.	Check the operation of the device.	
One or both ionization LEDs do <b>not</b> light up.	The ionization LEDs are malfunctioning.	Check the operation of the device.	8.2, 11

### 8.2 Checking the Operation of the Device

#### Procedure

- ▶ Switch the device off and on again.
- ▶ If the problem persists: The device is defective.
- ▶ Switch off the device.
- ▶ Disconnect the device from the power supply at the installation site.
- ▶ Secure the device to prevent it from restarting.
- ▶ Contact Sartorius Service.

## 9 Decommissioning

### 9.1 Decommission the Device

#### Procedure

- ▶ Switch off the device.
- ▶ Disconnect the device from the power supply.
- ▶ Clean the device.

## 10 Storage and Shipping

### 10.1 Decontaminating the Device

The device does **not** contain any hazardous materials that necessitate special disposal measures. If the device has come into contact with hazardous substances: Steps must be taken to ensure proper decontamination and declarations.

#### Procedure

- ▶ Decontaminate the device. The operator of the device is responsible for adhering to local government regulations on the proper decontamination and declaration for transport and disposal.

### 10.2 Storage

#### Procedure

- ▶ Clean the device.
- ▶ Store the device according to the ambient conditions (see chapter 12.2, page 13).

### 10.3 Returning the Device and Parts

Defective devices or parts can be sent back to Sartorius. Returned devices must be clean, decontaminated, and properly packed.

Transport damage as well as measures for subsequent cleaning and disinfection of the device or parts by Sartorius are charged to sender.

Devices contaminated with hazardous materials, e.g., harmful biological or chemical substances, will **not** be accepted for repair or disposal. The devices must be decontaminated before shipping (for decontamination, see chapter "10.1 Decontaminating the Device", page 12).

#### Procedure

- ▶ Decommission the device.
- ▶ Contact Sartorius Service for instructions on how to return devices or parts (please refer to [www.sartorius.com](http://www.sartorius.com)).
- ▶ Pack the device and its parts properly for return.

## 11 Disposal

### 11.1 Disposing of the Device and Parts

The device and the device accessories must be disposed of properly by disposal facilities.

The packaging consists of environmentally friendly materials that are intended to be used as secondary raw materials.

#### Procedure

- ▶ Dispose of the device. Follow the disposal instructions on our website ([www.sartorius.com](http://www.sartorius.com)).
- ▶ Dispose of the packaging in accordance with local government regulations.

## 12 Technical Data

### 12.1 Dimensions and Weights

	Unit	Value
Dimensions		
Length x width x height	mm	120 x 120 x 55
Weight	g	478

### 12.2 Ambient Conditions

	Unit	Value
Indoor use only		
Pollution level, according to EN 61010-1		2
Height above sea level, maximum	m	2000
Suitable for protection class		
Device protection class according to DIN EN 60529-1		IP40
Temperature		
In operation	°C	-20 - +40
During storage and transport	°C	-20 - +40
Relative air humidity		
<b>Non-condensing</b>	%	20 - 90

### 12.3 Electrical Specifications

#### 12.3.1 Power Supply

	Unit	Value
Power supply only permitted using Sartorius power supply unit		
Primary		
AC voltage	V	100 - 240
Frequency	Hz	50 - 60
Power Consumption	W	2

	Unit	Value
Secondary		
DC voltage	V	12
Current consumption, nominal	mA	250
Current consumption, maximum	mA	800
Power supply consumption	W	0.85
Overvoltage category according to IEC 60,664-1		II
For further data, see label on the power supply unit		
Ground fault circuit interrupter present at the installation site		
Socket at the installation site with fuse, maximum	A	16

### 12.3.2 Safety and Electromagnetic Compatibility

Safety of electrical equipment, according to EN 61326-1 Electrical equipment for measurement, control, and laboratory use – EMC requirements – Part 1: General requirements (IEC 61326-1)

Interference immunity: Suitable for use in industrial areas

Transient emissions

Class B: Suitable for use in residential areas and areas that are directly connected to a low voltage network that also supplies residential buildings.

### 12.3.3 Neutralizing

	Unit	Value
Neutralizing	V	± 30
Ion source: High-voltage ion generator		

## 12.4 Operating Area

	Unit	Value
Effective removal of objects (in front of the device)	mm	30 – 60

## 12.5 Discharging Times


Fig. 3: Discharging times (mm/sec) for neutralizing an electrostatically charged, insulated conductor | neutralizing from 1000 V to 100 V

## 12.6 Cleaning Agents and Cleaning Processes

---

Housing: Wipe with damp cloth with rubbing alcohol, maximum 70%

---

# 13 Sartorius Service

Sartorius Service is at your disposal for queries about the device. Please visit the Sartorius website ([www.sartorius.com](http://www.sartorius.com)) for information about the service addresses, services provided, or to contact a local representative.

When contacting Sartorius Service with questions about the system or in the event of malfunctions, be sure to have the device information, e.g., serial number, close at hand. This information can be found on the manufacturer's ID label.

# 14 Conformity Documents

The attached documents declare the conformity of the device with the designated directives or standards.

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Sartorius Lab Instruments GmbH & Co. KG  
г. Гёттинген, 13.12.2021

---

Dr. Reinhard Baumfalk  
Head of Product Development (LPS)

---

Halil Yildirim  
Product Compliance Officer (SLI)

Sartorius Lab Instruments GmbH & Co. KG  
Otto-Brenner-Strasse 20  
37079 Goettingen, Germany

Phone: +49 551 308 0  
[www.sartorius.com](http://www.sartorius.com)

The information and figures contained in these instructions correspond to the version date specified below.

Sartorius reserves the right to make changes to the technology, features, specifications and design of the equipment without notice.

Masculine or feminine forms are used to facilitate legibility in these instructions and always simultaneously denote all genders.

Copyright notice:

These instructions, including all components, are protected by copyright.

Any use beyond the limits of the copyright law is not permitted without our approval.

This applies in particular to reprinting, translation and editing irrespective of the type of media used.

Last updated:

12 | 2021