

# Virosart® Max\*

Maximize your virus  
filter performance with  
adsorptive pre-filtration



## Product Information

Virosart® Max is a specifically optimized virus pre-filter significantly increasing downstream virus filter performance. This filter combines size exclusion mechanism with efficient adsorptive capacities to increase the robustness of the following virus filter. As a result, Virosart® Max ensures the highest protection of your final virus retentive membrane, significantly increasing its robustness and capacity.

### Description

Choose your perfect fit from the Sartorius virus clearance strategy summarizing orthogonal technologies, manufacturing solutions, validation support and consultancy.

The orthogonal technologies from Sartorius consist of virus inactivation as well as virus removal by chromatography and virus filtration. The Virosart® product range includes four

different virus retentive membranes, in order to provide the best solution for every application. The performance of the final virus filter mainly depends on the feed stream properties. Therefore, appropriate feed stream conditioning i.e. through efficient pre-filtration is a major aspect within the development of a virus filtration step.

### **Application & Positioning**

The main application for Virosart® Max is the protection of the final virus retentive membrane. For that purpose this filter is used at the end of the purification process in-line with the final virus filtration step of the biopharmaceutical product.

The optimized pre-filter to final-filter ratio should be determined during development of the virus filtration step.

### **Product Benefits**

Virosart® Max provides highest protection of the virus retentive membrane with significant increase in robustness and capacity of the final virus filter. The most challenging molecules for the final virus retentive membrane are aggregates and/or small hydrophobic molecules. Virosart® Max binds aggregates very efficiently through hydrophobic interactions, independently of process conditions such as conductivity. The clean triple layer membrane material provides highest adsorptive capacities optimized for this process step.

With the T-Style design the Virosart® Max Maxicaps® is ideal for easy installation of multiple filters in series or parallel.

### **Integrity Testing**

Virosart® Max filters are tested for integrity using a water-based integrity test with the Sartocheck® technology of Sartorius Stedim Biotech.

\* The patented technology (DE 10 2011 105 525 B4) binds aggregates efficiently through hydrophobic interactions with polyamide, independently of process conditions such as conductivity from biological feed streams (mAbs, plasma derivatives or recombinant proteins).

# Technical Data




	Minisart® & Sartolab	Capsule & Midicaps®		T-Style Maxicaps® & Cartridge
Nominal filtration area	5 & 21 cm <sup>2</sup>	220 cm <sup>2</sup> 0.24 ft <sup>2</sup>	0.21 m <sup>2</sup> 2.3 ft <sup>2</sup>	0.7 m <sup>2</sup> , 1.4 m <sup>2</sup> , 2.1 m <sup>2</sup> 7.5 ft <sup>2</sup> , 15 ft <sup>2</sup> , 22.6 ft <sup>2</sup>
To be used for	<ul style="list-style-type: none"> <li>Scale-down work</li> <li>Filtration work for capacity studies</li> <li>Optimizing the pre-filter final-filter-ratio of the final virus filtration step</li> </ul>	<ul style="list-style-type: none"> <li>Scale-up studies</li> </ul>		<ul style="list-style-type: none"> <li>Large scale manufacturing</li> </ul>
Typical filtration volume	< 200 mL	< 80 L		> 80 L
Available connectors	<ul style="list-style-type: none"> <li>Female luer lock inlet &amp; male luer lock outlet</li> </ul>	<ul style="list-style-type: none"> <li>Midicaps® &amp; Capsule: ¾" triclamp (sanitary) connector inlet &amp; outlet</li> </ul>		<ul style="list-style-type: none"> <li>Maxicaps®: Sanitary inlet &amp; outlet adapter</li> <li>Cartridges: S-adapter top, 2 flange bayonet adapter with double o-ring bottom</li> </ul>
Sterilization	<ul style="list-style-type: none"> <li>Autoclaving: 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul>	<ul style="list-style-type: none"> <li>Autoclaving: 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul>		<ul style="list-style-type: none"> <li>Maxicaps®: Autoclaving 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul>
	<p>⚠ No inline steaming of Minisart®</p>	<p>⚠ No inline steaming of Capsule &amp; Midicaps®</p>		<p>⚠ No inline steaming of Maxicaps®</p> <ul style="list-style-type: none"> <li>Cartridges: Steaming 121°C @ 2.0 bar   29 psi for 30 min up to 2 cycles</li> </ul>
Operating parameters	<ul style="list-style-type: none"> <li>In the direction of filtration: max. 5.0 bar   72.5 psi at 20°C; max. 0.5 bar   7.3 psi at 121°C</li> <li>In the reversed direction of filtration: max. 0.5 bar   7.3 psi at 20°C</li> </ul>			
Water based diffusion test at 2.0 bar   29 psi	N/A	<ul style="list-style-type: none"> <li>2 mL/min (220 cm<sup>2</sup>)</li> <li>10 mL/min (0.27 m<sup>2</sup>)</li> </ul>		<ul style="list-style-type: none"> <li>8 mL/min (0.7 m<sup>2</sup>)</li> <li>16 mL/min (1.4 m<sup>2</sup>)</li> <li>24 mL/min (2.1 m<sup>2</sup>)</li> </ul>

# Materials

## Device

Cartridge, T-Style Maxicaps® , Midicaps® & Capsule	
Supportive Fleece	Polypropylene
Core (not capsule)	Polypropylene
End Caps	Polypropylene
Capsule Housing	Polypropylene
Minisart®	
Housing	Polypropylene

Sartolab	
Housing	Polycarbonate
Membrane	
Material	Optimized polyamide
Pore size	0.1 µm (nominal)
Format	Triple layer



Construction of Viroart® Max cartridge and capsule with zoom on cross section and membrane.


# Performance

## Characteristic Water Flow Rates

The flow rate achieved through the Virosart® Max, in-line with the final virus filter, is determined by the flow rate of the final virus filter!


Virosart® Max 30" Cartridges & 30" T-Style Maxicaps®

(2.1 m<sup>2</sup> | 22.6 ft<sup>2</sup>)




Virosart® Max 20" Cartridges & 20" T-Style Maxicaps®

(1.4 m<sup>2</sup> | 15 ft<sup>2</sup>)




Virosart® Max 10" Cartridges & 10" T-Style Maxicaps®

(0.7 m<sup>2</sup> | 7.5 ft<sup>2</sup>)




Virosart® Max Midicaps®

(0.21 m<sup>2</sup> | 2.3 ft<sup>2</sup>)



Virosart® Max Capsules

(220 cm<sup>2</sup> | 0.24 ft<sup>2</sup>)




Regulatory Compliance

- Each individual filter is tested for integrity (except 54AMI-----B and 54ASA-----V)
- Designed, developed and manufactured in accordance with an ISO 9001 certified Quality Management System
- Meet or exceed the requirements for WFI quality standards set by the current USP
- Non pyrogenic according to USP Bacterial Endotoxins
- USP Plastic Class Test VI

Technical References

Validation Guide	2650008-000-00
Extractables Guide	SPK5803-e

Ordering Information




**Cartridge**

54A	25	58	■
-----	----	----	---

**Adapter**  
25: S-adapter top,  
1 flange bayonet  
adapter with double  
O-ring bottom

**Filter size**  
N1: 10" 0.7 m<sup>2</sup> | 7.5 ft<sup>2</sup>  
N2: 20" 1.4 m<sup>2</sup> | 15 ft<sup>2</sup>  
N3: 30" 2.1 m<sup>2</sup> | 22.6 ft<sup>2</sup>




**T-Style Maxicaps®**

54A	83	58	■	--	SS
-----	----	----	---	----	----

**Filter size**  
N1: 10" 0.7 m<sup>2</sup> | 7.5 ft<sup>2</sup>  
N2: 20" 1.4 m<sup>2</sup> | 15 ft<sup>2</sup>  
N3: 30" 2.1 m<sup>2</sup> | 22.6 ft<sup>2</sup>

**Adapter**  
SS: 1 ½" Tri-Clamp (sanitary)  
inlet & outlet adapter



**Midicaps® & Capsules**


54A	■	58	■	--	FF	--	■
-----	---	----	---	----	----	----	---

13: Capsules  
53: Midicaps®

**Filter size**  
N4: 220 cm<sup>2</sup> | 0.24 ft<sup>2</sup>  
N9: 0.21 m<sup>2</sup> | 2.3 ft<sup>2</sup>

**Connector**  
FF: ¾" Tri-Clamp (sanitary)  
inlet & outlet adapter

**Units per package**  
B: Five pieces  
V: Two pieces



**Minisart®**

54A	■	--	--	--	--	--	■
-----	---	----	----	----	----	----	---

**Filter size**  
M1: Minisart®  
SA: Sartolab

**Units per package**  
B: Five pieces  
V: Two pieces

# Accessories & Services

## Virus Filtration

Virosart® HF is a high speed virus filter for mAbs and recombinant proteins. The filter combines highest virus safety with superior capacities.



## Integrity Testing using Sartocheck®

Fully automated Virosart® integrity testing to guarantee intactness of the Virosart® filter applying pre- and post-use diffusion tests.



## Ready-to use Filter Transfer Sets


Simplify your daily routine work by using modular filter assembly.

## Single-use Systems

Flexible processing with FlexAct® VR system for production from pilot plants up to commercial processing.

## Customized Systems

High level of automation and individual requirements can be realized by customized single-use or hybrid solutions.



Sartorius Confidence® Virus Clearance Services are the perfect complement to Virosart® virus filters. Our services provide:

- Virus clearance studies
- Process design support
- Optimization support

We use a variety of different relevant and model viruses including MVM, MuLV, Reo-3 and HSV-1. The combination of product and services provides you with a comprehensive virus clearance solution that gives you the confidence you need to proceed.

## BioOutsource Testing Services

Your partner to assure virus safety for your process by MCB | WCB characterization, bulk harvest testing.

**Germany**

Sartorius Stedim Biotech GmbH  
August-Spindler-Strasse 11  
37079 Goettingen  
Phone +49 551 308 0

**USA**

Sartorius Stedim North America Inc.  
565 Johnson Avenue  
Bohemia, NY 11716  
Toll-Free +1 800 368 7178



For further contacts, visit  
[www.sartorius.com](http://www.sartorius.com)